



Development of quantum sensor technology for synchrotron applications

Lothar Weinhardt, Mary Blankenship, Kilian Bitsch, Constantin Wansorra, Ralph Steininger, Thomas Spangenberg, Dirk Hauschild, Nik Arldt, Tonya Vitova, Sebastian Kempf, and Clemens Heske

Institute for Photon Science and Synchrotron Radiation (IPS), Karlsruhe Institute of Technology (KIT)

Institute for Chemical Technology and Polymer Chemistry (ITCP), Karlsruhe Institute of Technology (KIT)




Department of Chemistry and Biochemistry, University of Nevada, Las Vegas (UNLV)




Electronic and chemical structure of applied material systems

Electron spectroscopies

- Charge-carrier transport
- Chemical bonding
- Interaction with light




X-ray spectroscopies




Weinhardt et al., *Advanced Materials* **31**, 1806660 (2019)


What can we learn with x-ray spectroscopy?



Soft x-ray emission spectroscopy (XES) is sensitive to the local chemical environment – example: sulfur



- XES probes the local partial density of states
“Valence band structure from the viewpoint of a selected core level (atom)”
- Extremely sensitive to local bonding environment








Phys. Stat. Sol. **187**, 13 (2001) (Cu₂S and CuS)
Appl. Phys. Lett. **83**, 2360 (2003) (HgS and ZnS)
J. Chem. Phys. **124**, 074705 (2006) (MoS₂)
Phys. Rev. B **75**, 165207 (2007), *Phys. Rev. B* **79**, 165305 (2009) (CdS)
Appl. Phys. Lett. **97**, 074101 (2010) (In₂S₃)
J. Phys. Chem. C **119**, 10412 (2015) (In₂S₃)
ACS Appl. Mater. Interf. **8**, 33256 (2016) (In₂S₃ and ZnS)
Anal. Chem. **23**, 8300 (2021) (CdSO₄)
Sol. RRL **7**, 2201091 (2023) (ZnS)
Prog. Photovolt.: Res. Appl. **32**, 904 (2024) (In₂S₃ and ZnS)
Phys. Chem. Chem. Phys. **26**, 26389 (2024) (ZnS, CdS, HgS, In₂S₃)
Adv. Mater. Interf. **12**, 2401002 (2025) (CdS, CdSO₄)

...

Cation speciation in sulfates



— Experiment

— Calculation





Weinhardt et al.,
 Analytical Chemistry **93**, 8300 (2021)

Resonant Inelastic X-ray Scattering (RIXS) map of C₆₀




- Before SALS: Measure spectra at a few select energies
- With SALS: Closely spaced absorption scan with a spectrum at every excitation energy




Weinhardt et al., JCP **135**, 104705 (2011)

2D/3D perovskite heterostructures for thin-film solar cells





- Detailed chemical picture
- Organic materials sensitive to x-ray irradiation
- dose must be kept low
- here by continuous moving sample under beam
→ exposure of a given point: 50 ms



Hauschild et al.,
ACS Appl. Mater. Interfaces **13**, 53202 (2021)

In situ / operando experiments

Gas or liquid continuously pumped through probing volume
Separated from UHV by a thin membrane (Si_3N_4 , SiC, C)



Fuchs et al., *PRL* **100**, 249802 (2008)







SALSA (2008)

Fuchs et al., *Nucl. Instrum. Methods Phys. Res.* **585**, 172 (2008)
Blum et al., *Rev. Sci. Instr.* **80**, 123102 (2009)
Weinhardt et al., *J. Electr. Spectr. Rel. Phenom.* **188**, 111 (2013)
Benkert et al., *Rev. Sci. Instr.* **85**, 015119 (2014)
Weinhardt et al., *J. Synchrotron Rad.* **28**, 609 (2021)

Operando study of catalytical processes

- Direct study of reactant, intermediate, and product species


Seitz et al., *J. Phys. Chem. C* **126**, 20998 (2022)



Fouda et al., *J. Phys. Chem. Lett.* **11**, 7476 (2020)


High efficiency soft x-ray spectrometer

- 1 Beamline flux
- 2 Solid angle
- 3 Efficiency of optics
- 4 Detector efficiency
- 5 Dead times





1

Beamline flux



X-SPEC @ KIT Light Source




UHV and *in situ* / *operando* environments
XES, XAS, RIXS, EXAFS, PES, HAXPES


Weinhardt et al., J. Synchrotron Rad. **28**, 609 (2021)

Dead times

X-SPEC @ KIT Light Source



- Coordinated continuous movement of monochromator and undulator
- Fast CMOS sensor

Weinhardt et al., J. Synchrotron Rad. **28**, 609 (2021)



Rapid RIXS maps (rRIXS) at X-SPEC (KIT Light Source)

1540 s

118 s



60 s

30 s



Weinhardt et al. *J. Synchrotron Rad.* **31**, 1481 (2024)


Quantum x-ray sensors



Krantz et al., Appl. Phys. Lett. 124, 032601 (2024)

High efficiency soft x-ray spectrometer


- 1 Beamline flux
- 2 Solid angle
- 3 Efficiency of optics
- 4 Detector efficiency
- 5 Dead times



High efficiency soft x-ray spectrometer


	Grating spectrometer (X-SPEC @ KIT Light Source)	Quantum sensor
1 Beamline flux		
2 Solid angle	0.0017 – 0.0045 % of ½ sphere	Much larger!
3 Efficiency of optics	2 – 25 %	??
4 Detector efficiency	> 90 %	> 90 %
5 Dead times		

Is this the best we can do?



Quantum x-ray sensors

Poster Mary Blankenship



Work items

- Develop suitable filters to block IR radiation but transmit (soft) x-rays
- We need many pixels for high count rates!
- Filter “unwanted” spectral components?


Ultra-thin Filter Foil Candidates

K. Bitsch et al., Master Thesis (2024)


Commercially available ultra-thin foils for filtering IR radiation:

- 50 nm Pd, 70 nm Al, 100 nm Zr
- 21.5 nm Al/47.6 nm polyimide

→ Fragility of these foils makes measurements challenging!



unbroken foil example




IR Radiation Transmission

Cryostat base temperatures reach as low as 7 mK, sensor array cooled to temperatures of 17 mK

M. Krantz et al., Appl. Phys. Lett. 124, 032601 (2024)

Detector resolution is noticeably reduced even at such low transmission

- Multiple filters are necessary to block enough of the IR radiation



K. Bitsch et al., Master Thesis (2024)

Summary

- Soft X-ray Spectroscopies are very powerful techniques for electronic and chemical structure studies
- Material research needs high-transmission x-ray spectrometers
- Quantum detectors promise the next big leap forward

